



**Giles Scientific Inc.**  
P.O. Box 4306  
Santa Barbara, CA 93140 USA  
Phone: (805) 963-3876  
Fax: (805) 963-7768  
Email: [giles@biomic.com](mailto:giles@biomic.com)  
[www.biomic.com](http://www.biomic.com)

## Press Release - January 2007

### 2007 CLSI Guidelines in BIOMIC System

The latest 2007 CLSI M100 & M-44A (bacteria/yeast antibiotic disk test) guidelines are included in the BIOMIC software. All test results are monitored by an advanced Expert System to flag unlikely results. Users may also define their own expert rules.

BIOMIC V3 Microbiology System also reads commercial bacteria and yeast identification (ID) panels including RapID™ (Remel), API® (bioMerieux), Crystal™ (BD), and CHROMagar™ (BD). High resolution color test images (both AST and ID) can be saved for further evaluation.

BIOMIC V3 serves as a primary or complimentary AST-ID system for microbiology laboratories of all sizes and setups. A bi-directional LIS interface can be added to further speed workflow. BIOMIC V3 connects to a standard Windows-based PC and requires no routine maintenance or service.

#### About Giles Scientific, Inc

Giles Scientific is a privately held company that has served the clinical microbiology and pharmaceutical laboratory markets since 1984. The Company develops and manufactures a line of digital imaging instrumentation based on proprietary technologies covered by numerous patents, copyrights and trademarks, and distributes these products in 50+ countries.

**Contact Giles Scientific Inc:** website [www.biomic.com](http://www.biomic.com) email: [giles@biomic.com](mailto:giles@biomic.com) or call 800-603-9290 for more information.