



Giles Scientific Inc.
P.O. Box 4306
Santa Barbara, CA 93140 USA
Phone: (805) 963-3876
Fax: (805) 963-7768
Email: giles@biomic.com
www.biomic.com

Press Release – May 2009

CAP Checklist added to BIOMIC V3

The BIOMIC V3 Microbiology System now contains a built-in CAP (College of American Pathologists) Bacteriology Checklist to facilitate laboratory inspections. This includes all bacteriology-specific CAP Questions with a direct link to its BIOMIC V3 test and report feature, including QC Summary Reports and Antibiotic Disk Lots Log. This new addition to the BIOMIC V3 Quality Control Center includes electronic documentation of all QC test results.

BIOMIC V3's color digital imaging technology automates the reading, interpretation, and reporting of CLSI, EUCAST, and BSAC antibiotic disk diffusion tests. The system also reads all varieties of API[®] (bioMerieux) RapID[™] (Remel) and Crystal[™] (BD) ID panels and automates colony counting. BIOMIC V3 serves as a primary or secondary AST-ID system for microbiology laboratories of all sizes and setups.

About Giles Scientific Inc:

Giles Scientific has served clinical, pharmaceutical and veterinary microbiology laboratories since 1984. The Company develops, manufactures, sells and supports our customers with our color digital imaging equipment in 55+ countries.

Contact Giles Scientific Inc:

Website: www.biomic.com

Email: giles@biomic.com

Phone: 805-963-3876