



**Giles Scientific Inc.**  
P.O. Box 4306  
Santa Barbara, CA 93140 USA  
Phone: (805) 963-3876  
Fax: (805) 963-7768  
Email: [giles@biomic.com](mailto:giles@biomic.com)  
[www.biomic.com](http://www.biomic.com)

## **Press Release – May 2008**

### **3M™ Petrifilm™ Automated Reading and Image Storage**

DigiCounter™ software now automatically counts all varieties of 3M™ Petrifilm™. Entire films or selected areas are instantly counted. Counts may be differentiated by size and or color, and touching colonies can be discriminated. Individual colonies may be manually included or excluded. High resolution color digital images are saved and can be reviewed any time.

DigiCounter software is included with both BIOMIC® V3 and TRINITY™ V3 microbiology systems from Giles Scientific (California, USA).

#### **About Giles Scientific, Inc**

Giles Scientific has served clinical, pharmaceutical and veterinary microbiology laboratories since 1984. The Company develops, manufactures, sells and supports our customers with our color digital imaging equipment in 50+ countries.

**Contact Giles Scientific Inc:** website [www.biomic.com](http://www.biomic.com) email: [giles@biomic.com](mailto:giles@biomic.com) or call 800-603-9290 to request information.