



**Giles Scientific Inc.**  
P.O. Box 4306  
Santa Barbara, CA 93140 USA  
Phone: (805) 963-3876  
Fax: (805) 963-7768  
Email: [giles@biomic.com](mailto:giles@biomic.com)  
[www.biomic.com](http://www.biomic.com)

## **Press Release – January 2009**

### **Multiple Language Translations in BIOMIC V3**

BIOMIC software is now available in Chinese, Spanish, French, and Portuguese. The addition of other languages is based on demand from our global network of BIOMIC V3 customers.

BIOMIC V3's advanced color digital imaging technology automates the reading, interpretation, and reporting of CLSI and BSAC antibiotic disk diffusion tests and commercial identification panels including RapID™ (Remel), API® (bioMerieux), BBL™ Crystal™ (BD). The addition of a bar code reader coupled with bi-directional LIS/LIMS interface capability further speeds workflow. BIOMIC V3 serves as a primary or secondary AST-ID system for microbiology laboratories of all sizes and setups.

#### **About Giles Scientific Inc:**

Giles Scientific has served clinical, pharmaceutical and veterinary microbiology laboratories since 1984. The Company develops, manufactures, sells and supports our customers with our color digital imaging equipment in 55+ countries.

#### **Contact Giles Scientific Inc:**

Website: [www.biomic.com](http://www.biomic.com)

Email: [giles@biomic.com](mailto:giles@biomic.com)

Phone: 805-963-3876