



Giles Scientific Inc.
P.O. Box 4306
Santa Barbara, CA 93140 USA
Phone: (805) 963-3876
Fax: (805) 963-7768
Email: giles@biomic.com
www.biomic.com

Press Release - June 2007

New Etest® Click and Read

The BIOMIC V3 Microbiology System now performs semi-automated reading for all varieties of Etest® (AB BIODISK) strips. An enlarged plate image with the Etest strip(s) and its drug specific MIC scale is displayed on-screen. The user selects the appropriate MIC value, triggering an S-I-R interpretation. Etest result(s) can be printed or transferred to an LIS. High resolution color Etest images can also be saved for future review.

The BIOMIC V3's advanced color digital imaging technology automates the reading, interpretation, and recording of CLSI / BSAC antibiotic disk diffusion tests and commercial identification panels including RapID™ (Remel), API® (bioMerieux), BBL™ Crystal™ (BD). The antibiotic resistance detection and quality control features follow current CLSI guidelines.

About Giles Scientific, Inc

Giles Scientific has served veterinary, clinical, and pharmaceutical microbiology laboratories since 1984. The Company develops, manufactures, sells and supports our customers with our color digital imaging equipment in 50+ countries.

Contact Giles Scientific Inc: website www.biomic.com email: giles@biomic.com or call 800-603-9290 to request information.