



Giles Scientific Inc.
P.O. Box 4306
Santa Barbara, CA 93140 USA
Phone: (805) 963-3876
Fax: (805) 963-7768
Email: giles@biomic.com
www.biomic.com

Press Release – November 2009

Microbial Water Testing with BIOMIC V3

BIOMIC V3 provides automated microbiological water testing including species identification, organism colony counting, and antibiogram phenotyping.

BIOMIC V3 identifies bacteria and yeast by reading and interpreting commercial identification panels (including all varieties of API[®] from bioMérieux). Colonies are counted on agar plates and can simultaneously identify species on chromogenic agar plates with visible and ultraviolet light. BIOMIC V3 also provides antibiogram phenotypes by reading and interpreting antibiotic disk diffusion zones, with 20 antibiotics markers used in one published study.

About Giles Scientific Inc:

Giles Scientific has served clinical, environmental, pharmaceutical and veterinary microbiology laboratories since 1984. The Company develops, manufactures, sells and supports our customers with our color digital imaging equipment in 55+ countries.

Contact Giles Scientific Inc:

Website: www.biomic.com

Email: giles@biomic.com

Phone: 805-963-3876