



**Giles Scientific Inc.**  
P.O. Box 4306  
Santa Barbara, CA 93140 USA  
Phone: (805) 963-3876  
Fax: (805) 963-7768  
Email: [giles@biomic.com](mailto:giles@biomic.com)  
[www.biomic.com](http://www.biomic.com)

## **Press Release – August 2009**

### **Automated Reading of Rosco Neo-Sensitabs™ with BIOMIC V3**

The BIOMIC V3 Microbiology System now automates the reading and interpretation of Rosco Neo-Sensitabs. Neo-Sensitabs comprise 125+ different antimicrobial agents. BIOMIC V3 interprets Neo-Sensitabs based on current CLSI and EUCAST guidelines.

BIOMIC V3's color digital imaging technology also reads and interprets all brands of disks for antibiotic susceptibility testing. Additional features include reading and interpreting of identification panels including API® (bioMérieux), RapID™ (Remel), Crystal™ (BD), chromogenic agar and colony count test plates. Bar code reading and a bi-directional LIS/LIMS interface further speeds testing workflow. BIOMIC V3 serves as a primary or secondary AST-ID system for microbiology laboratories.

#### **About Giles Scientific Inc:**

Giles Scientific has served clinical, pharmaceutical and veterinary microbiology laboratories since 1984. The Company develops, manufactures, sells and supports our customers with our color digital imaging equipment in 55+ countries.

#### **Contact Giles Scientific Inc:**

Website: [www.biomic.com](http://www.biomic.com)

Email: [giles@biomic.com](mailto:giles@biomic.com)

Phone: 805-963-3876