



Giles Scientific Inc.
P.O. Box 4306
Santa Barbara, CA 93140 USA
Phone: (805) 963-3876
Fax: (805) 963-7768
Email: giles@biomic.com
www.biomic.com

Press Release – June 2019

Bottom Camera Reader accessory for BIOMIC® V3

Giles Scientific now offers a BIOMIC V3 optional accessory providing automated 96-well microtiter plate* reading from the bottom of the plate, versus reading from the top with BIOMIC V3. The Bottom Camera Reader utilizes the same software and monitor as BIOMIC V3. Reading 96-well microtiter plates from the bottom provides more accurate test results by avoiding interference including condensation and bubbles on the top of the wells. All plate images are saved in high resolution and available for future review.

About Giles Scientific Inc:

Giles Scientific has served clinical and industrial microbiology laboratories since 1984. Giles develops, manufactures, sells and supports our customers in 100 countries.

Contact Giles Scientific Inc:

Website: www.biomic.com

Phone: 1-805-963-3876

*96-Well Microtiter - Note in USA: For Research Use Only. Not for use in diagnostic procedures.



Giles Scientific Inc.
P.O. Box 4306
Santa Barbara, CA 93140 USA
Phone: (805) 963-3876
Fax: (805) 963-7768
Email: giles@biomic.com
www.biomic.com

Press Release – June 2019

Bottom Camera Reader accessory for BIOMIC® V3

Giles Scientific now offers a BIOMIC V3 optional accessory providing automated 96-well microtiter plate* reading from the bottom of the plate, versus reading from the top with BIOMIC V3. The Bottom Camera Reader utilizes the same software and monitor as BIOMIC V3. Reading 96-well microtiter plates from the bottom provides more accurate test results by avoiding interference including condensation and bubbles on the top of the wells. All plate images are saved in high resolution and available for future review.

About Giles Scientific Inc:

Giles Scientific has served clinical and industrial microbiology laboratories since 1984. Giles develops, manufactures, sells and supports our customers in 100 countries.

Contact Giles Scientific Inc:

Website: www.biomic.com

Phone: 1-805-963-3876

*96-Well Microtiter - Note in USA: For Research Use Only. Not for use in diagnostic procedures.