

TRINITY™ V3

Industrial Microbiology System

Automate and standardize zone reading and colony counting

- Instant digital reading & calculations
- IQ/OQ/PQ validation documents
- 24 inch touch screen
- Live customer support 7 days per week
- No routine maintenance required



Optional Testing Modules:



Antibiotic Potency Zone Reading


- USP, EP, JP methods
- 21 CFR Part 11 compliant



Colony Counting

- Colony Count Plates
- Chromogenic Agar Plates

Optional TRINITY™ V3 Accessories:



Flat & Beveled Peni Cylinders



Peni Cylinder Dispenser