

Press Release – April 2024

2024 CLSI and EUCAST Guidelines in the BIOMIC® V3 Microbiology System

The following guidelines are updated in the 2024 BIOMIC V3 software:

CLSI Guidelines:

- M100, 34th Edition: Performance Standards for Antimicrobial Susceptibility Testing (February 2024)
- VET01S: Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated from Animals, 7th Edition (January 2024)

EUCAST Guidelines:

- Breakpoint Tables for Interpretation of MICs and Zone Diameters v 14.0 (January 2024)
- Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. v 14.0 (January 2024)
- Zone diameter breakpoints for rapid antimicrobial susceptibility testing (RAST) directly from blood culture bottles. v 6.1 (June 2023)

About Giles Scientific Inc:

Giles Scientific (Santa Barbara, California, USA) has served clinical, environmental, pharmaceutical and veterinary microbiology laboratories since 1984. Giles Scientific develops, manufactures, sells and supports customers with our digital imaging equipment in 100+ countries. Visit www.biomic.com to learn more.