

BIOMIC® V3

Clinical Microbiology System

Open system to automate and standardize microbiology tests

- Instant digital reading & interpretation
- 2-way LIS/LIMS interface
- CLSI and EUCAST guidelines
- Save high resolution test images
- Resistance monitoring & antibiograms
- 24 inch touch screen monitor
- No routine maintenance



Optional Testing Modules:



Disk Diffusion

- BD, Oxoid, etc
- Etest® & MIC Test Strips



96-Well Microtiter*

- Custom Plates
- Commercial Plates



Colony Counting

- Colony Count Plates
- Chromogenic Agar Plates



Organism Identification

- API® (bioMerieux)
- RapID™ (Thermo)
- Crystal™ (BD)
- Liofilchem® ID Systems



MIC Strip

- Etest® (bioMerieux)
- MTS® (Liofilchem)

*96-Well Microtiter - Note in USA: For Research Use Only. Not for use in diagnostic procedures.



Quality Control Testing & Inventory Management

- AST, ID & reagent QC testing
- QC summary reports
- Documentation of disks, plates, reagents & media