

BIOMIC® V3

Clinical Microbiology System

Open system to automate and standardize microbiology tests

- Instant digital reading & interpretation
- 2-way LIS/LIMS interface
- CLSI and EUCAST guidelines
- Save high resolution test images
- Resistance monitoring & antibiograms
- 24 inch touch screen monitor
- No routine maintenance



Optional Testing Modules:




Disk Diffusion

- BD, Oxoid, etc




Organism Identification

- API® (bioMerieux)
- RapID™ (Thermo)
- Liofilchem® ID Systems




Agar Dilution*

- Standard Methods
- Custom Methods



Broth Microdilution*

- Custom Plates
- Commercial Plates



Colony Count

- Colony Count Plates
- Chromogenic Agar Plates



MIC Strip

- Etest® (bioMerieux)
- MTS® (Liofilchem)



Urine Screen

- Custom Methods
- Commercial Methods

*Agar Dilution & Broth Microdilution - Note in USA: For Research Use Only. Not for use in diagnostic procedures.