



# BIOMIC<sup>®</sup> V<sup>3</sup>

## Microbiology System for Veterinary Labs

Open system to automate and standardize manual microbiology tests

- Instant digital reading & interpretation
- 2-way LIS/LIMS interface
- CLSI veterinary guidelines
- Customized veterinary drug panels
- Resistance monitoring & antibiograms




### Testing Modules:



#### Antibiotic Susceptibility

- Disk Diffusion Tests
- Etest<sup>®</sup> & MIC Test Strips



#### 96-Well Microtiter\*

- Custom Plates
- Commercial Plates



#### Organism Identification


- API<sup>®</sup> panels
- RapID<sup>™</sup> panels
- Crystal<sup>™</sup> panels



#### Colony Counting

- Colony Count Plates
- Chromogenic Agar Plates

\*96-Well Microtiter - Note in USA: For Research Use Only. Not for use in diagnostic procedures.



### Quality Control Testing & Inventory Management

- AST, ID & reagent QC testing
- QC summary reports
- Documentation of disks, plates, reagents & media