

## Press Release – October 2025

### **Automated Code Recognition for Liofilchem® Antibiotic Disks with BIOMIC® V3**

BIOMIC V3 now automatically recognizes antibiotic disk codes from Liofilchem (Italy) antibiotic disks, in addition to BD and Oxoid disks. For example, Liofilchem disk code AM-10 is automatically recognized as Ampicillin. If a BIOMIC V3 customer uses disks from a different manufacturer (not Liofilchem, BD, Oxoid) the customer would simply assign the disk positions on the plate during a one-time setup.

#### **About Giles Scientific Inc:**

Giles Scientific (Santa Barbara, California, USA) has served clinical, environmental, pharmaceutical and veterinary microbiology laboratories since 1984. Giles Scientific develops, manufactures, sells and supports customers with our digital imaging equipment in 100+ countries. Visit [www.biomic.com](http://www.biomic.com) to learn more.