





BIOMIC[®] V³

Quality Control & Inventory Management


Instant Digital Reading & QC Testing




Antibiotic Disks



Etest[®]/MIC Test Strips



ID Panels



96-Well Microtiter Plates*

*96-Well Microtiter - Note in USA: For Research Use Only. Not for use in diagnostic procedures.

Inventory Management

Record:

- Entry Date
- Material
- Quantity
- Manufacturer
- Lot Number
- Expiration Date
- In-Use Date
- Out-Use Date
- Received By
- Visual Inspection
- Comments
- Custom Defined Fields

Date Received	Material	Manufacturer	Lot #	Expire Date	Comment	Last QC
6/11/2012	Amikacin	Remel	95865879909855	4/25/2013		8/22/2012
6/11/2012	Amoxicillin-clavulanate	Remel	857475785755443	3/27/2013		8/6/2012
6/11/2012	Ampicillin	bioMérieux	954788897677765	3/21/2013		8/6/2012
1/9/2012	Aztreonam	BD	789796416351591	2/27/2013		8/22/2012
1/9/2012	Cefepime	Remel	348754765453504	3/20/2013		8/22/2012
1/9/2012	Ceftazidime	Remel	15764354325543	5/20/13		8/22/2012
6/11/2012	Cephalexin	BD	784564534678344	3/13/2014		8/6/2012
1/9/2012	Chloramphenicol	BD	87857864784786	3/13/2013		8/22/2012
1/9/2012	Ciprofloxacin	Remel	987851654155488	7/31/2013		8/22/2012
1/9/2012	Colistin	Remel	985467545475454	1/2/2012		8/22/2012
6/11/2012	Concomitant	BD	96548498989656	6/12/2012		8/6/2012
1/9/2012	Meropenem	Remel	564675485698765	3/13/2013		8/22/2012
6/11/2012	Norfloxacin	Remel	46545436868758	2/12/2014		8/6/2012
1/9/2012	Piperacillin-tazobactam	BD	454895785455676	8/28/2013		8/22/2012
1/9/2012	Ticarcillin-clavulanate	BD	897654684165455	4/18/2014		8/22/2012
1/9/2012	Tobramycin	Remel	748854757654887	2/13/2014		8/22/2012

QC Summary Reports



- QC Cumulative Summary
- QC Detailed Summary
- QC Abbreviated
- QC Zone Trend Graph
- QC Comment Summary
- QC Annual Trend Graph
- MIC QC Cumulative Summary
- New Lot QC Summary
- ID Panel QC Summary
- New Lot ID Panel QC Detailed Summary
- New Lot ID Panel QC Summary
- Reagent QC Summary
- Thymidine Control Test Summary