



## 96-Well Microtiter\* (MP-MOD)

### Step 1:

#### Insert plate & select 'Read Microplate'

- Sensititre (TREK) AST\* & ID Plates
- MicroScan (Siemens) AST\* & ID Plates
- Custom 96-Well AST Plates\* (all brands)




21 inch  
Touch Screen


\*96-Well Microtiter - Note in USA: For Research Use Only. Not for use in diagnostic procedures.

### Step 2:


#### Review plate image & well readings



TREK YeastOne Y-09\*



Custom 96-Well\*



TREK RAPMYCO\*

### Step 3:

#### Transfer results to LIS/LIMS

	Drug Name	Raw	SIR	MIC	Expert
1	MIC Cefixor	S		0.50	
2	MIC Fluconazole	I		4.0	
3	MIC Tylosin	N/A		> 20	
4	MIC Sulphadimethoxime	R		> 256	
5	MIC Erythromycin	R		> 4.0	
6	MIC Gentamicin	S		<=0.50	
7	MIC Clindamycin	N/A		> 4.0	
8	MIC Entrofloxacin	S		<=0.13	
9	MIC Oxleytacycline	N/A		2.0	
10	MIC Neomycin	N/A		<=2.0	
11	MIC Sulphafurazole	I/I			
12	MIC Novobiocin	I/I			



Antibiograms

