


Peni Cylinder Dispenser & Peni Cylinders

Accessories to **TRINITY™ V3**



Peni Cylinder Dispenser

- Holds 300+ cylinders
- Dispense 4 or 6 cylinders per 90-100mm Petri plates
- Used with plastic or glass Petri plates
- Made of stainless steel & aluminum
- Dimensions: 13 x 13 x 31 cm



Flat Face or Beveled

Peni Cylinders

- Made of stainless steel
- Dimensions:
Inner diameter: 6mm
Outer diameter: 8mm
Length: 10mm
- Tolerance: ± 0.1 mm