



Giles Scientific Inc.
P.O. Box 4306
Santa Barbara, CA 93140 USA
Phone: (805) 963-3876
Fax: (805) 963-7768
Email: giles@biomic.com
www.biomic.com

Press Release – October 2013

Remote Access for BIOMIC® V3 Microbiology System

Customers can now remotely access BIOMIC V3 from additional computers in the lab network. This new feature provides significant labor savings as customers can instantly review test and QC results, reports, and plate images from any office or workstation in the lab.

BIOMIC V3 completes the automation in the microbiology lab by reading antibiotic disk diffusion following CLSI/EUCAST/BSAC guidelines, Etest (bioMerieux), MIC Test Strips (Liofilchem) and 96-well microtiter plates. BIOMIC V3 also automates API® (bioMerieux), RapID™ (Remel), BBL™ Crystal™ (BD), Liofilchem® ID Systems, colony count and chromogenic agar plates. Bar code entry and a 2-Way LIS/LIMS interface optimizes lab workflow. BIOMIC V3 serves as a primary or secondary AST-ID system for microbiology laboratories.

About Giles Scientific Inc:

Giles Scientific has served clinical, environmental, pharmaceutical and animal health microbiology laboratories since 1984. The Company develops, manufactures, sells and supports our customers in 75 countries.

Contact Giles Scientific Inc:

Website: www.biomic.com

Email: giles@biomic.com

Phone: 1-805-963-3876