

# TRINITY™ V3

## Antibiotic Potency Zone Reading

Automate and standardize zone reading and calculations for antibiotic potency assays

- Instant digital reading & calculations
- IQ/OQ/PQ validation documents
- 21 CFR Part 11 compliant
- 24 inch touch screen
- Live customer support 7 days per week




### Methods:

- USP 81 (United States Pharmacopeia)
- EP (European Pharmacopeia)
- BP (British Pharmacopeia)
- JP (Japanese Pharmacopeia)
- Large Plate

Zone Measurements (mm)						
Assay Title	Reference			Standard		
User Name	Zone	1	2	3	4	5
A	1	14.0	15.7	15.0	15.2	11.0
	2	14.0	12.4	15.0	12.2	15.1
	3	14.0	12.0	15.0	12.1	15.0
B	1	14.0	13.0	14.8	13.5	14.5
	2	14.0	12.8	14.9	14.0	15.7
	3	14.0	12.8	14.9	14.0	15.7
C	1	14.0	14.4	14.5	14.8	14.4
	2	14.0	14.0	14.9	14.8	14.4
	3	14.0	14.2	14.3	14.5	14.7
D	1	14.0	17.4	15.1	17.6	14.7
	2	14.0	17.4	14.8	17.4	15.2
	3	14.0	17.1	14.5	17.9	15.1
Unknown	1	14.0	13.0	14.8	13.5	14.5
	2	14.0	12.8	14.9	14.0	15.7
	3	14.0	12.8	14.9	14.0	15.7
2	1	14.0	14.4	14.6	14.8	14.4
	2	14.0	14.5	14.5	14.9	14.4
	3	14.0	14.2	14.3	14.6	14.7

### Optional TRINITY™ V3 Accessories:



Flat Face & Beveled Peni Cylinders



Peni Cylinder Dispenser