

Press Release – January 2025

New Model Upgrades for BIOMIC® V3 and TRINITY™ V3

The new hardware models are now available for the BIOMIC V3 and TRINITY V3 Microbiology Systems. A summary of updates is listed below:

- 20 megapixel camera replaces previous 5 megapixel camera
- 70% more linear pixel resolution
- 20% lower distortion levels
- 24 inch high resolution touch screen
- Windows 11 compatible
- New streamlined calibration procedure

About Giles Scientific Inc:

Giles Scientific (Santa Barbara, California, USA) has served clinical, environmental, pharmaceutical and veterinary microbiology laboratories since 1984. Giles Scientific develops, manufactures, sells and supports customers with our digital imaging equipment in 100+ countries. Visit www.biomic.com to learn more.